

February 28th, 2022
Rev 5

EOVE

Products Catalogue



Contents

| | |
|---------------------------------------|-----------|
| EOVE-150 | 2 |
| Packs | 3 |
| Bags | 6 |
| Accessories | 7 |
| List of references | 10 |
| EOVE-70 | 15 |
| Pack | 16 |
| Accessories | 16 |
| List of references | 17 |
| Common references | 18 |
| Accessories | 19 |
| List of references | 20 |
| Services | 22 |
| ANNEX : Order Form Instruction | 23 |

EOVE-150

Life Support Ventilator



PACKS

Life support ventilator EOVE-150

EO-150VNT EUR
EO-150VNT US

Thanks to its design and its high definition intuitive graphic interface, the EOVE-150 ventilator combines comfort and simplicity of use for caregivers to high ventilation performances. The vast choice of ventilation modes and circuit configurations ensure an effective ventilation for adult or pediatric patients from 3.5 kg.

The life support ventilator EOVE-150 pack includes :

- The Life Support ventilator EOVE-150
- A double limb adaptor
- A leakage / mouthpiece single limb adaptor
- An O2 inlet connector
- An AC power charger module (EUR or US according to the geographic area)
- A transport bag



EOVE 150



EO-DB2-1P-KIT



EO-LMPADAPT



EO-O2CON



EO-PWRCHRG



EO-CARBAG1X0

External battery pack

EO-BATPCK



The external battery pack offers 9 supplemental ventilation hours for an enhanced mobility and safety. A USB slot also allows the patient to connect an external tablet or a Smartphone in order to keep a graphic interface in all circumstances.

The external battery pack includes :

- An external battery
- A discharge cable to link the external battery to the EOVE-150
- An AC power charger for the external battery



EO-BAT9



EO-BATCHRG



EO-BATCBL

Transport bag for EOVE-150

EO-CARBAG1X0



The Transport bag is used to store and transport the EOVE-150 ventilator. Its design allows the patient to store the ventilator and its accessories, power source and circuit.

Travel bag for EOVE-150

EO-TRVELBAG1X0-EVO



The Travel bag allows the patient to move with the complete EOVE-150 ventilator. Two external battery bags (optional) can be hooked, giving the patient up to 18 hours of battery life.

Nomad bag for EOVE-150


EO-NOMADBAG-EVO



The Nomad bag allows to use only the ventilation module while moving, allowing the patient to move easily with his ventilator over the shoulder or as a backpack, keeping a graphic interface thanks to the Smartphone application or to a dedicated external tablet.

Battery bag for external battery pack

EO-BATBAG



Combined to the Nomad or Travel bag, the battery bag allows the patient to easily carry one or two external batteries, offering a mobility solution with 9 or 18 hours battery life.

Battery packs coupling cable

EO-CPLPACK



Thanks to two external battery coupling cable, the patient will get 18 hours of external battery life.

Coupling cable for AC power

EO-YCBLPWR



The Y cable AC power charger allows to link the AC power charger to the EOVE-150 external battery for more safety. In case of AC power loss, the ventilator will run on external battery, giving the patient more autonomy.

Car lighter DC cable

EO-CARCBL



The car lighter DC cable allows the patient to connect the EOVE-150 ventilator or a stand-alone ventilation module to a car lighter, for highly secure long distance travels.

Warning : for countries outside Europe

Car DC/AC box

EO-CARBOX



The car DC/AC box allows the patient to connect the EOVE-150 ventilator or a stand-alone ventilation module to a car lighter, for highly secure long distance travels.

Warning : for countries inside Europe

O2CEL CBL

FiO₂ kit cable



The FiO₂ kit cable allows to link the EOVE-150 ventilator or a stand-alone ventilation module to a FiO₂ measuring cell (not provided), and to activate associated alarms and record FiO₂ data.

Proximal flow sensor

EO-PFLOWS



Put at the end of the circuit and linked to the EOVE-150 ventilator or a stand-alone ventilation module, the proximal flow sensor allows to measure the expired volumes in single limb configuration with exhalation valve. Used in pediatric mode, it allows an effective and secure ventilation for patients from 3.5 kg, and allows the volume monitoring from 30 ml.

Remote alarm cable (2 m) for Monnal alarm box

EO-RMOTCBL2



The 2 meters free ended nurse call cable allows to link the EO-150 ventilator to nurse call at hospital and to warn caregivers in case of alarm, so the patient's safety can be guaranteed.

Remote alarm cable (15 m) for Monnal alarm box

EO-RMOTCBL15



The 15 meters free ended nurse call cable allows to link the EO-150 ventilator to nurse call at hospital and to warn caregivers in case of alarm, so the patient's safety can be guaranteed.

Kit double limb adaptor

EO-DB2-1P-KIT



It contains the double Limb adaptor and the interface. It could be used to switch from a EO-DBADAPT to the new Double Limb Adaptor

List of reference

| PACKS | REFERENCE | PHOTO |
|--|-------------------------------|---|
| <i>Life Support ventilator EOVE-150</i> | EO-150VNT US EO-150VNT EUR |  |
| <i>External battery pack</i> | EO-BATPCK |  |
| ACCESSORIES | REFERENCE | PHOTO |
| <i>Transport bag for EOVE-150 ventilator</i> | EO-CARBAG1X0 |  |
| <i>Travel bag for EOVE-150 ventilator</i> | EO-TRVELBAG1X0-EVO |  |
| <i>Nomad bag for EOVE-150 ventilation module</i> | EO-NOMADBAG-EVO |  |
| <i>Battery bag for external battery</i> | EO-BATBAG |  |
| <i>Supplemental belt</i> | EO-BAGBELT |  |

| | | |
|--|--|---|
| <i>Double branch adaptor V1</i> | EO-DBADAPT |  |
| <i>Double Branch Adaptor V2 Single Patient</i> | EO-DB2-1P |  |
| <i>Kit Evolution Double Branch Adaptor V2 Single Patient</i> | EO-DB2-1P-KIT |  |
| <i>Double Branch Adaptor V2 Single Patient Module Interface</i> | EO-DB2-1P-INTERFACE |  |
| <i>Exhalation valve kit (membrane + cover) for double limb adaptor</i> | EO-EXHVALVE |  |
| <i>Screw for double limb adaptor</i> | DBSCREEW |  |
| <i>Double limb adaptor protection</i> | EO-DBPROTEC |  |
| <i>Leakage / Mouthpiece adaptor</i> | EO-LMPADAPT EO-LMPADAPT5 (x5 packaging) |  |
| <i>O₂ inlet connector</i> | EO-O2CON |  |

| | | |
|---|--------------|---|
| <i>External battery</i> | EO-BAT9 |  |
| <i>External battery discharge cable</i> | EO-BATCBL |  |
| <i>External battery charger</i> | EO-BATCHRG |  |
| <i>Coupling cable for AC power / external battery</i> | EO-YCBLPWR |  |
| <i>Car lighter DC cable</i> | EO-CARCBL |  |
| <i>FiO₂ kit cable</i> | O2CELCBL |  |
| <i>Proximal flow sensor</i> | EO-PFLOWS |  |
| <i>15 m remote alarm cable for Monnal alarm box</i> | EO-RMOTCBL15 |  |
| <i>2 m remote alarm cable for Monnal alarm box</i> | EO-RMOTCBL2 |  |

| EOVE-150 COMPONENTS | REFERENCE | PHOTO |
|--|------------------|--|
| <i>Docking station with EO-Display</i> | EO-DCK1SLT |  |
| <i>Ventilation module</i> | EO-VM150 |  |
| <i>Air inlet standard filter (x10)</i> | EO-AIFF |  |
| <i>Air inlet pollen filter (x10)</i> | EO-AIPF |  |

EOVE-70

Secretion management device



PACK

Secretion management device EOVE-70

EO-70SMD EUR
EO-70SMD US

The Pack "Secretion management device EOVE-70" includes :

- The secretion management device EOVE-70
- A transport bag
- A charger module



ACCESSORIES

Foot switch EOVE-70

EO-70FSWITCH



The foot switch helps caregivers to control the inhale et exhale cycle with a foot to keep the hands available to make recruitment maneuver on the patient.

LIST OF REFERENCES

PACK**REFERENCE****PHOTO**

*Secretion management device
EOVE-70*

EO-70SMDEUR



ACCESSORIES**REFERENCE****PHOTO**

EO-70 Foot switch

EO-70FSWITCH



Common references to EOVE-150 & EOVE-70



ACCESSORIES

Nonin XPod for SpO2 monitoring

EO-SP02CBL



The Nonin Xpod SpO2 cable allows to link the EOVE-150 ventilator or the EOVE-70 secretion management device to a Nonin SpO2 sensor, and to activate associated alarms and record SpO2 data.

Trolley for EOVE devices

EO-TROLLEY



The EOVE Trolley allows to put EOVE device in a roll stand at the hospital and to move it easily.

LIST OF REFERENCES

| ACCESSORIES | REFERENCE | PHOTO |
|--|------------|---|
| <i>Nonin XPod for SpO2 monitoring</i> | EO-SPO2CBL |  |
| <i>Upright brackets</i> | EO-UPRIGHT |  |
| <i>Trolley</i> | EO-TROLLEY |  |
| <i>AC power charger module for EOVE-150 / EOVE-70</i> | EO-PWRCHRG |  |
| <i>Docking / Module locking screw for EOVE-150 / EOVE-70</i> | DCKSCREW |  |

